20000017440

HB eco_092018

IKA

HB eco

Betriebsanleitung Ursprungssprache	DE	4
Operating instructions	EN	10
Mode d'emploi	FR	16
Instrucciones de uso	ES	22
Istruzioni per l'uso	IT	28
Руководство пользователя	RU	34
Manual de instruções	PT	40
Instrukcja eksploatacji	PL	46
ليغش تالام لعت	AR	52
Kullanma talimatları	TR	58

IKA-Werke, Germany
Reg. No. 004343

Geräteaufbau/Display - Device setup/Display

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
Fig. 1

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Contents

	Page		
Device setup/Display	2	Interfaces and outputs	13
Declaration of conformity	10	Maintenance and cleaning	13
Explication of warning symbols	10	Error codes	14
Safety instructions	11	Warranty	15
Correct use	11	Technical Data	15
Unpacking	12		
Commissioning	12		

Declaration of conformity

We declare under our sole responsibility that this product corresponds to the directives 2014/35/EU, 2014/30/EU and 2011/65/EU and conforms with the following standards or normative documents: EN 61010-1, EN 61010-2-010, EN 61326-1, EN 60529 and EN ISO 12100.

Explication of warning symbols

Indicates an (extremely) hazardous situation, which, if not avoided, will result in death, serious injury.

Indicates a potentially hazardous situation, which, if not avoided, can result in death, serious injury.

Indicates a potentially hazardous situation, which, if not avoided, can result in injury.

Indicates practices which, if not avoided, can result in equipment damage.

DANGER - note on hazards arising from a hot surface.

Safety instructions

For your protection

- Read the operating instructions completely before starting up and follow the safety instructions.
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the device.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- Wear your personal protective equipment in accordance with the hazard category of the medium to be processed. Otherwise there is a risk of splashing liquids.
- Set up the device in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- Prior to each use, always check the device for damage. Do not use damaged components.

CAUTION

Only process and heat up media that has a flash point higher than the adjusted safe temperature limit of the heating bath that has been set.

DANGER

Risk of burns! During operation, the heating bath housing can get hot.

- When emptying the device use only the handles to carry and hold it.
- Prior to filling or emptying the heating bath, the device must be switched off and disconnected from the power supply at the plug.
- Only fill or empty the heating bath when it is cold.
- Empty the heating bath prior to transporting it.
- Never operate the heating bath without tempering medium.

DANGER

The preferred tempering medium in the heating bath is water. There is a risk of burning when using tempering media with lower flash-points!

- Before use, calculate the optimum filling level of the tempering medium! Pay special attention to the change in volume caused by heating and the displacement that occurs when immersing an object, for example, an evaporating flask.

- When using the heating bath in combination with a rotary evaporator, the heating bath temperature must not be allowed to rise to a value higher than the boiling point of the solvent at normal pressure, since if the evaporating flask glass were to break there would be a hazard due to liquid spraying out (for instance breakage of the evaporating flask glass during distillation of ether when the medium temperature is at 60 °C).

CAUTION

When working with the IKA rotary evaporator, be aware of a hazard arising from breakage of the evaporating flask glass.

CAUTION

Be aware of a hazard due to lack of grip on a wet evaporating flask.

- The use of demineralised water is recommended.
- Beware of hazards due to flammable materials.
- Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
- Do not operate the device in explosive atmospheres, with hazardous substances or under water.
- The device will automatically restart in mode C following any interruption to the power supply.
- The device can only be disconnected from the mains supply by pulling out the mains plug or the connector plug.
- The socket for the mains cord must be easily accessible.

For protection of the equipment

- The voltage stated on the type plate must correspond to the mains voltage.
- Socket must be earthed (protective ground contact).
- Protect the device and accessories from bumps and impacts.
- The device may only be opened by experts.

Solvents

- Solvents can be hazardous to health. Therefore comply with the relevant warnings and refer to the relevant safety data sheet (Internet).

Correct use

• Use

The IKA tempering baths HB eco is laboratory device and is suitable for directly tempering substances filled into the bath container. They are also suitable for indirect tempering of substances filled in glass containers when the glass containers are immersed in the actual tempering medium. A rotary glass container is particularly advantageous, for example, when used in conjunction with an IKA rotary evaporator.

CAUTION

The device is not intended to prepare food!

• Area of use (only indoors)

- Laboratories
- Pharmacies
- Schools
- Universities

This device is suitable for use in all areas except:

- Residential areas
- Areas that are connected directly to a low-voltage supply network that also supplies residential areas.

The safety of the user cannot be guaranteed:

- If the device is operated with accessories that are not supplied or recommended by the manufacturer;
- If the device is operated improperly or contrary to the manufacturer's specifications;
- If the device or the printed circuit board are modified by the third parties.

Unpacking

• Unpacking


- Unpack the device carefully,
- Any damage should be notified immediately to the shipping agent (post, rail or logistics company).

• Scope of delivery

- Heating bath **IKA** HB eco
- Power supply cable
- USB cable
- Operating instructions
- Warranty card

Commissioning

Observe the ambient conditions (temperature, humidity, etc.) listed under "Technical data".

The unit is ready for service when the mains plug has been plugged in.

Prior to commissioning the device, determine the filling quantity of the tempering medium due to the used evaporator piston size (approx. 2.5 l tempering medium by using a standard 1 l evaporator piston).

- Switch on the device using the mains switch (B).
- At every start, the display will show all display segments (Fig. 2), the software version (Fig. 3), and the operating mode (Fig. 4).
- The heating bath temperature is held constant by the control

circuit and is additionally monitored by the safety circuit.

In the event of an error occurred in the control circuit, the heating bath is switched off permanently by the safety circuit.

An error in the control or safety circuits is shown on the display. The heating function can no longer be started.

- When switching off the heating function and the medium temperature is higher than 50 °C, the display will alternately show the set temperature (Fig. 5) and "HOT" (Fig. 9). Otherwise, the display will alternately show the set temperature (Fig. 5) and OFF (Fig. 8).

When switching on the heating function, the display will alternately show set temperature (Fig. 5) and the actual temperature (Fig. 6 or Fig. 7).

Setting the target temperature

Activate the SET function by turning and pressing the rotating/pressing knob (C).

- Set the desired heating bath temperature by turning the rotating/pressing knob (C).

- The LED setpoint temperature (H) will be always on whenever the set temperature is shown. (Fig. 5)
- Confirm the value by pressing the rotating/pressing knob (C).
- Start the heating process by pressing the rotating/pressing knob (C).

Setting the operating modes

Operating mode A

- Switch on the device at the mains switch (B).
- The heating function is switched off.
- The target value is set to 20 °C.
- The target value can be adjusted.
- The heating function must be restarted following a power failure.
- Mode A appears on the display when switching on the device.

Factory setting: Mode A

Operating mode B

- Switch on the device at the mains switch (B).
- The heating function is switched off.
- The target value is set to 20 °C or the most recent set temperature.
- The target value can be adjusted.
- The heating function must be restarted following a power failure.
- Mode B appears on the display when switching on the device.

Operating mode C

- Switch on the device at the mains switch (B).
- The heating function is switched on or off depending on the last selected setting.
- The last target value set in mode B applies.
- The target value cannot be adjusted.
- The heating function will be automatically restarted following a power failure.
- Mode C appears on the display when switching on the device.

Switching the operating modes

- Switch the device off at the mains switch (B).
- Press and hold rotating/pressing knob (C) and switch the device on at the mains switch (B). Release the rotating/pressing knob (C) after 2 seconds.
- Sequence A, B, C, A etc.
- The operating mode appears on the display when switching on the device.

Safety temperature

The safety circuit prevents the temperature of the heating bath from getting too high as a result of controller error. Once the safety temperature has been reached, the device switches off permanently.

Furthermore, dry running of the heating bath is detected. The device switches off permanently.

The dry running function detects both accidental heating of the heating bath without liquid and dry running caused by water evaporation starting at a set temperature of 60 °C. The error message E 26 (Fig. 11) appears and the heating bath is switched off permanently. To remedy this error, see "Error codes".

Controlling the medium temperature

The medium temperature is limited by the set safety temperature controlled by means of a PID controller and measured by the PT 1000 temperature sensor, and the medium is heated up as quickly as possible without overshooting to the set temperature. The PID controller adapts to the different tempering media and ensures optimum temperature control with minimum temperature drift and fluctuation.

Optimum control is only achieved when the medium is mixed using a rotating evaporator piston.

- Set the desired medium temperature between room temperature and 99 °C by turning the rotating/pressing knob (C).
- Start the heating function by pressing the rotating/pressing knob (C).
- The LED heating (G) appears.
- The heating bath is heated up to the set temperature.
- The set and actual temperature referred to the medium are shown on the display.

Interfaces and outputs

Note: Please comply with the system requirements together with the operating instructions and help section included with the software.

USB interface

The Universal Serial Bus (USB) is a serial bus for connecting the device to the PC. Equipped with USB devices can be connected to a PC during operation (hot plugging). Connected devices and their properties are automatically recognized. Use the USB interface to update the firmware.

USB device drivers

First, download the latest driver for **IKA** devices with USB interface from:

<http://www.ika.com/ika/lws/download/usb-driver.zip>

Install the driver by running the setup file. Then connect the **IKA** device through the USB data cable to the PC. The data communication is via a virtual COM port.

USB cable A - B


This cable is used to connect the USB interface to a PC.

Maintenance and cleaning

The device is maintenance-free. It is only subject to the natural wear and tear of components and their statistical failure rate.

Cleaning

Remove the device from the mains before cleaning.

Only use cleaning materials recommended by **IKA**:

Dirt	Cleaning agent
Dyes	Isopropyl alcohol
Building materials	Water containing detergent/isopropyl alcohol
Cosmetics	Water containing detergent/isopropyl alcohol
Food	Water containing detergent
Fuels	Water containing detergent
Other materials	Please consult IKA

Wear protective gloves during cleaning the devices.

Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.

Do not allow moisture to get into the device when cleaning.

Before using another than the recommended method for cleaning or decontamination, you must ascertain with **IKA** that this method does not destroy the device.

Ordering spare parts

When ordering spare parts, please give:

- Device type,
- Serial number, see rating plate,
- Position number and description of spare part, see www.ika.com,
- Software version.

Repairs

Please only send devices in for repair that have been cleaned and are free of materials which might present health hazards.

For this, use the “**certificate of compliance**” form which you can obtain from **IKA** or can download a version for printing from the **IKA** website at www.ika.com.

If your appliance requires repair, return it in its original packaging. Storage packaging is not sufficient when sending the device - also use appropriate transport packaging.

Error codes

The fault is shown by an error code on the display as following if the error occurs.

Proceed as follows in such cases:

- Switch the device off with mains switch,
- Carry out corrective measures,
- Restart the device.

Error code	Causes	Effect	Solutions
E 3	- Temperature inside the device is too high	Heating off	- Switch off the device - Let the device cool down - Switch on the device
E 9	- Fault in storing the safety circuit set temperatures - Memory chip (EPROM) defective	Heating off	- Switch off the device - Let the device cool down - Switch on the device
E 21	- Safety relay does not open	Heating off	- Switch off the device and switch on again
E 24	- Set safety temperature exceeded	Heating off	- Switch off the device - Let the device cool down - Switch on the device
E 25	- The heater control circuit switch (TRIAC) is faulty. The heater or the supply line is disconnected.	Heating off	- Switch off the device and switch on again
E 26	- Dry running	Heating off	- Switch off the device - Let the device cool down - Fill up the tempering medium - Switch on the device
E 27	- Error in calibration	Heating off	- Switch off the device and switch on again
E 28	- Breakage of the controller sensor	Heating off	- Switch off the device and switch on again
E 29	- Short-circuit in the safety sensor	Heating off	- Switch off the device and switch on again
E 30	- Short-circuit in the controller sensor	Heating off	- Switch off the device and switch on again
E 31	- Breakage of the safety sensor	Heating off	- Switch off the device and switch on again
E 32	- Temperature deviation is too big	Heating off	- Switch off the device and switch on again

If the actions described fails to resolve the fault or another error code is displayed then take one of the following steps:

- Contact the service department,
- Send the device for repair, including a short description of the fault.

Warranty

In accordance with **IKA** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine directly to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

Technical data

Heating output	W	1350
Heating temperature range	°C	Room temperature...99
Heat control		LED
Set temperature resolution	K	±1
Filling volume max.	l	4
Filling point min.	mm	60
Material in contact with medium		Stainless steel 1.4404
Fixed safety temperature	°C	109
Safety class DIN 12876		I
Outer height	mm	190
Inner height	mm	130
Dimensions (W x H x D)	mm	330 x 190 x 325
Weight	kg	3.9
Permissible ambient temperature	°C	5...40
Permissible relative humidity	%	80
Protection class according to DIN EN 60529		IP 21
RS 232 interface		no
USB interface		yes
Analog output		no
Voltage	V	200...240 ±10 % 100...120 ±10 %
Frequency	Hz	50/60
Power input	W	1350
Power input standby	W	3

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IKA-Werke GmbH & Co.KG