


Catalogue No.	Cone Size	Dia "D" <sub>≈</sub>	L1±2	L2±2	L3±2	Stopcock Bore Size
GPS04820202	14/23	9	60	45	45	2mm
GPS04820202X	14/23	9	60	45	45	2.5mm



## GLASS APPS

<b>GLASS/APPS</b>		MATERIAL: BOROSILICATE GLASS, 3.3 EXPANSION, TYPE I, HGA 1	
SCALE : NTS	NAME	SIGN	DATE
DRN.	DRS		
CHKD.	ANIL VERMA		
TITLE: Adapter, Flushing 14/23 With T-Bore PTFE Stopcock, Cone		PROJ.:	DRAWING No.
		QUANTITY	REV.
		SHEET	0

OLD REF. DRG. NO. ASSEMBLY DRG. NO.	SR.NO. DATE	REVISION	DRN. APD.	1	2	3	4	5	6
--	----------------	----------	--------------	---	---	---	---	---	---

ALL DIMENSIONS ARE IN MM, UNLESS OTHERWISE SPECIFIED