


Catalogue No.	Cone Size	Dia 'D'≈	L1±2	L2±2	L3±2	Stopcock Bore Size
GPS04820201	24/40	9	75	45	45	2mm
GPS04820201X	24/40	9	75	45	45	2.5mm



SCALE : NTS

MATERIAL:
BOROSILICATE GLASS,
3.3 EXPANSION, TYPE I,
HGA 1

GLASS APPS

TITLE: Adapter, Flushing 24/40

With T-Bore PTFE Stopcock, Cone

DRN.	NAME	SIGN	DATE
CHK.D.	ANIL VERMA		

DRN.	DRS	PROJ.:
CHK.D.	ANIL VERMA	

DRN.	DATE	QUANTITY	SHEET	DRAWING No.	REV.
CHK.D.				GPS0482020XX-1	0

OLD REF. DRG. NO.					
ASSEMBLY DRG. NO.					
SRNO.	DATE	REVISION	DRN.	APPD.	

ALL DIMENSIONS ARE IN MM, UNLESS OTHERWISE SPECIFIED